

Comparison of binding and functional data on a set of safety-related GPCRs



Jolas, T¹, Migeon, J², Fouchaq, B¹

¹Eurofins Cerep SA, Le bois l'Evêque, 86600 Celle l'Evescault, France

²Eurofins Panlabs Inc, 15318 NE 95th Street, Redmond WA 98052, USA

Eurofins Pharma Discovery Services

Introduction

Getting a better idea of the safety fingerprint of lead compounds early in the discovery process may reduce the risk of off-target interactions leading to adverse drug reactions (ADRs) (Bowes *et al.*, 2012). *In vitro* pharmacology profiling is a classical approach that has been used for this purpose.

For a set of 282 to 616 pharmacologically relevant compounds, we compared binding and functional data obtained on several GPCRs selected for their cardiovascular ADRs. Those GPCRs were part of the targets listed by Bowes *et al.* The objective of this study was to compare the sensitivity of the two different methods using dose-response determinations and assess how the binding and functional assays were related to one another.

Methods

Target	# cpds	Function measure	G protein	Agonist	Radioligand
5-HT _{1A}	465	cAMP	Gi	8-OH-DPAT	8-OH-DPAT (agonist)
5-HT _{2B}	616	Ca ²⁺ flux	Gq	serotonin	DOI (agonist)
α _{1A} Adrenoceptor	494	Ca ²⁺ flux	Gq	epinephrine	prazosin
α _{2A} Adrenoceptor	395	Impedance	Gi	epinephrine	RX821002
α _{2B} Adrenoceptor	490	cAMP	Gi	dexmedetomidine	RX821002
A _{2A} Adenosine	371	cAMP	Gs	NECA	CGS21680 (agonist)
β ₁ Adrenoceptor	282	cAMP	Gs	isoproterenol	CGP12177 (agonist)
β ₂ Adrenoceptor	424	cAMP	Gs	isoproterenol	CGP12177 (agonist)
D ₁ Dopamine	540	cAMP	Gs	dopamine	SCH 23390
H ₁ Histamine	460	Ca ²⁺ flux	Gq	Histamine	pyrilamine
H ₂ Histamine	491	cAMP	Gs	Histamine	APT
M ₂ Muscarinic	517	cAMP	Gi	Acetylcholine	AF-DX384

Radioligand binding was assessed at equilibrium, using the ligands indicated in the table, at recombinantly expressed human receptors.

Functional assessment was by:

- HTRF - *homogeneous time resolved fluorescence*
- Ca²⁺ flux – *fluorescent calcium sensitive dye*
- CDS – *cellular dielectric spectroscopy, impedance*

In functional cell-based assays, compounds were assessed for effect on their own (compound agonism) and for ability to inhibit the response to the reference agonists shown in the table (compound antagonism). Compounds were tested over 8 concentrations from a maximum of 100µM or 30µM.

Compounds were selected using the BioPrint database (Eurofins Cerep) as being marketed or withdrawn drugs and pharmacological reference compounds

Results

Radioligand Binding Assays are more sensitive than Cell-Based Assays More antagonists are detected than agonists




Figure 1. Number of concentration-response curves generated using Radioligand Binding Assays (RBA) (IC₅₀) and the functional Cell-Based Assays (CBA) (IC₅₀ for antagonists and EC₅₀ for agonists) across a panel of 12 GPCR targets

Abbreviations:
5HT_{1A}, serotonin 5-HT_{1A} receptor
5HT_{2B}, serotonin 5-HT_{2A} receptor
α_{1A}, α_{1A} adrenoceptor
α_{2A}, α_{2A} adrenoceptor
α_{2B}, adenosine A_{2A} receptor
A_{2A}, adenosine A_{2B} receptor
β₁, β₁ adrenoceptor
β₂, β₂ adrenoceptor
D₁, dopamine D₁ receptor
H₁, histamine H₁ receptor
H₂, histamine H₂ receptor
M₂, muscarinic M₂ receptor

- For each target tested, radioligand binding assays gave the greatest number of actives (defined IC₅₀ value).
- In functional assays, generally, more antagonists than agonists were detected; 5-HT_{1A} was the single exception.

Agonists → Binding IC50s and functional EC50s do not correlate for most targets




Figure 2. Comparison of pIC₅₀ values from RBA (x-axis) and pEC₅₀ from CBA (y-axis) for a subset of BioPrint compounds across a panel of 12 GPCR targets.

Compounds inactive in one or other assay are shown with an arbitrary pIC₅₀ or pEC₅₀ value of 1 (blue dots).

- Correlation between binding IC₅₀s and functional EC₅₀s with a R²>0.6 was obtained in 2 out of the 12 targets

Antagonists → Binding IC50s and functional IC50s correlate for most targets




Figure 3. Comparison of pIC₅₀ values from RBA (x-axis) and pIC₅₀ from CBA (y-axis) for a subset of BioPrint compounds across a panel of 12 GPCR targets.

Compounds inactive in one or other assay are shown with an arbitrary pIC₅₀ value of 1 (blue dots).

- Correlation between binding and functional IC₅₀s with a R²>0.6 was obtained in 8 out of the 12 targets.
- For targets with good correlation, between 88% and 100% of the test compounds with a binding IC₅₀ below 1µM have a calculated IC₅₀ in the functional assay.

Antagonists → Screening at 10 µM does not detect as many active compounds




Figure 4. Comparison of binding (x-axis) and antagonist properties (y-axis) for a subset of BioPrint compounds across a panel of 12 GPCR targets. Compounds were tested at 10 µM.

Black lines indicate 25% inhibition that in average is considered to be the limit of noise level.

Numbers in the upper corner of each graph indicate number of compounds with binding and functional effects above 25% / number of compounds with binding and functional effects above 50% versus number of compounds with binding and functional IC₅₀s

- Single concentration screening does not detect as many compounds active in binding and inhibition functional assays when compared to testing at several concentrations.
- Comparison between binding and functional data is more complex when performing single concentration screening.

Conclusions

- Radioligand binding assays detect more interactions between compound and receptor than functional cell-based assays.
- For most targets, a poor correlation was found between agonist EC₅₀ and binding IC₅₀ whereas a good correlation was found between antagonist IC₅₀ and binding IC₅₀, but correlation is a property of the individual tests and the functional mode being run.
- Most test compounds (88-100%) showing strong binding properties (IC₅₀ < 1µM) have a calculable IC₅₀ in the functional CBA assay. For binding IC₅₀s higher than 1 µM, a population of compounds had binding interactions without detectable cellular effects. Similarly, a population of compounds had agonist or antagonist effects without binding interactions. Those populations need to be further assessed for selectivity of the response, sensitivity of the assay and allosteric modulation.

In conclusion,

- affinity (binding) and functional (agonist/antagonist) data,
- multi-concentration testing,

are required to provide an integrated assessment of compound interaction with secondary targets.

References

Bowes *et al* (2012) Reducing safety-related drug attrition: the use of *in vitro* pharmacological profiling. *Nat Rev Drug Discov.* 2012 Dec;11(12):909-22

Acknowledgments

We wish to thank the Secondary Pharmacology group at Astra Zeneca for their contribution to this work and Jenny Thouvenin for poster preparation