



FASTEST TURNAROUND TIME FOR LESS MONEY!

Today, synthetic genes have become affordable and easy to order.

The highly automated work flow of Eurofins MWG Operon allows the fastest turnaround times on the market.

Service Overview

- Fastest delivery times
- Great price
- Unbeatable quality
- Excellent support
- Convenient online services

We offer standard
gene synthesis

at 0.36 €/bp in just 8 days*

without any additional
express fee expenses.

Save time and money with gene synthesis and get your perfect gene without any cloning on your part!

Service Overview & Easy and Convenient Online Ordering

Eurofins MWG Operon is an expert in gene synthesis for many years and is able to provide you any clone and sequence optimised for gene expression in any organism you need.

Service Overview

- Gene synthesis service - standard and complex genes in any size
- Gene Evolution libraries – customised design
- Site-directed mutagenesis – single or multiple mutations, deletions and insertions
- Plasmid preparation – from 1µg up to 10 mg plasmid DNA

Easy and Convenient Online Ordering

Our gene synthesis order menu is easy to use. After entry, sequences are analysed by our GENEius online tool regarding GC content, repetitive structures and hairpin structures. GENEius is also able to optimise the codon usage for any type of organism.

Sequences applicable for our standardised synthesis processes can then be ordered directly. If submitted sequences turn out to be potentially complex the sequences will be analysed manually by our experts.

You can also use the order tool to request a quote first before you decide to order. All orders, quotes and the tracking of orders can be archived and monitored in your online account.

*Standard genes up to 1000 bp, moderate GC content (40-65%), without extensive homopolymer stretches (> 20 bp), without extensive repeats (> 20 bp) or hairpin structures.

Local currencies: 0.33 GBP/bp -2.94 NOK/bp- 3.66 SEK/bp – 2.70 DKK/bp - 0.53 CHF/bp

¹ GENEius software designed and developed for Eurofins MWG Operon by BioLink Informationstechnologie GmbH.

Any Questions?

Feel free to contact our Customer Support team if you need any assistance.



Email:
support@eurofinsdna.com

Phone:
+49 8092 8289-77

Our official business hours are:
8 a.m. – 6 p.m. CET

Free Phone Numbers:

Italy	0800 785950
France	0800 903807
Denmark	8088 1262
UK	0800 0323135
Ireland	1800 555056
Sweden	020 798148
Norway	0800 13844
Finland	0800 112744
Switzerland	0800 562013