



Global Quality Control

Global and competitive Quality control (QC)

Eurofins offers one of the largest networks of quality control laboratories operating in the USA, Germany, France, Switzerland, Italy, Slovakia, Spain and Denmark. The presence of our local subsidiaries in the main markets is essential to deliver our outstanding performance with low turn around times.

Complete support for your product development cycle

We have the largest battery of different analytical and microbiological methods developed for quality control of raw materials, finished products and water control. Our main areas of expertise are:

- API, impurities testing
- Stability studies
- packaging testing (suitability, protection, sterility)
- cleaning validation
- environmental documentation
- sterility testing
- hygiene monitoring

Our broad experience in these different fields allows us to tailor our services to the customer requirements and to support the product development from the beginning up to marketing.



Eurofins is a life sciences company operating laboratories that provide a comprehensive range of analytical testing services, drawing on the latest developments in biotechnology. The Group specialises in delivering analytical testing and advisory services to clients from a wide range of industries including the pharmaceutical, food and environmental sectors.