



Instrumentation – Denmark

Liquid Chromatography

Dionex Ultimate HPLC	2
Thermo Ultimate 3000RS UPLC	2
HPLC Detectors	
Dionex Refractive Index Detector Shodes RI-101 (RI)	1
Dionex Photodiode Array Detector UVD340U (PDA)	1
Dionex VWD-3100 Detector	1
Dionex Diode Array Detectors DAD-3000 (DAD)	3
Thermo VWD-3100 Detector	1
Thermo Photoluminescence Detector FLD-3000 (FLD)	1

Gas Chromatography

Agilent 6890A GC-FID	2
Perkin Elmer Headspace Sampler TurboMatrix 40 HS	2

Ion Chromatography

Thermo Integrion	1
Dionex IC ICS-1600	1

Metal Analysis

Thermo ICP-AES ICAP 6300 and 7400	2
Milestone Microwave Oven Ethos	1
Milestone Ultrawave	1
Anton Paar Microwave Oven Multitwaves	2
Agilent ICP MS 7700x	1
Agilent ICP-MS 7900	3

Spectrophotometry

Shimadzu UV-1800	1
Thermo FT-IR Spectrofotometer Nicolet IS10	1

Microbiology

Termarks Incubators	15
Memmert Incubators	7
Laf Benches	8
Geringe Autoclave	1
Biotek Instruments Microplate Readers	2
Cold Storage Rooms (5°C)	2
AB 3500 XL Genetic Analyzer (MicroSeq)	1
Veriti Thermal Cyclers	2
Binder Incubator	1
Lab refrigerators	8
Lab Freezers, -80°C	2
Lab Freezers, -20°C	2


Microbiology (cont.)

Labguard for temperature and humidity monitoring	1
Maldi-TOF Microflex LT/SH	1
Spectrophotometer	1

Stability

Sanyo PSC062.XHX.C	3
Binder KBF 720, 25°C/60%RH	2

Miscellaneous

Shimadzu TOC-V CPH	1
Shimadzu TOC-V CSH	1
Methrohm Autotitrator 702 SM Titrino	1
Anton Paar Densitometer DMA38	1
B&S Polarimeter ADP220	1
Mettler Seven Compact S220	1
Mettler Seven Multi	1
Mettler Toledo Karl Fisher Titrator DL32 Coulometric	1

Metrology

Mettler Toledo AT 261 Delta range	1
Mettler XP205	3
Mettler XP603S DR	1

For equipment details, see VMP.